

X Implant CUBE

Basic Implant Engine

- 9 programmable memories
- Automatic overload protection function
- Actual RPM display describing working condition
- Motor auto-calibration function
- Ergonomic foot control pedal



ACL(B) - 41I / 42I / 43I / 45I



CRB26XX / CRB27XX



FORTE100αEI

Specification

Contra Angle Handpiece

	CRB26XX	CRB27XX	ACL(B)-41I	ACL(B)-42I	ACL(B)-43I	ACL(B)-45I
Gear Ratio	20 : 1	32 : 1	20 : 1	16 : 1	64 : 1	32 : 1
Max. Torque	55 N·cm	55 N·cm	55 N·cm	55 N·cm	65 N·cm	65 N·cm
Max. RPM	2,500	1,560	2,500	3,100	780	1,560
Min. RPM	30	20	30	40	10	20
Weight	84g	84g	60g	60g	64g	60g

Control Box

Input	AC 100-120 V, 220-240 V
Output	230 W
Fuse	5 A / 2 A
Frequency	50/60 Hz
Maximum Pump	Max .75mℓ/min
Dimensions (mm)	W 205 x D 210 x H 136
Weight	3.3kg

Foot Controller

